

APBA REGION 6 SPRING 2021 SPRING MEETING

FEBRUARY 20/21, 2021

**PLACE: HOLIDAY INN
27355 CARRONADE DRIVE
PERRYSBURG, OHIO 43551**

PHONE: 1-(567) 331-8900

RATES:\$99 plus tax

DIRECTIONS: Exit 193 I-75; east on US 20; hotel on the right across from the old HIFQ.

NOTES:

1. Please inform the Reservation Desk that you are attending the APBA Region 6 Meeting. There are 20 rooms being held for Friday night; 27 rooms for Saturday night. The hotel will hold our rooms until **February 5, 2021**.
2. The Region 6 Inboard Safety Team (Box 21) will have the "Test Cell" at poolside from 7:00am to 10:30am per a sign-up sheet schedule to be posted on Friday afternoon. Walk-ins cannot be given testing due to restrictions in place to meet COVID-19 requirements. **A \$20 donation is greatly appreciated!**
3. There will be no awards and presentations this year.

AGENDA – SATURDAY, FEBRUARY 20, 2021

7:00 a.m.	SAFETY CELL (until 9:00 a.m. or later if needed)	POOLSIDE
1:00 p.m.	OFFICIATING MEETING	CARRONADE ROOM 2
3:30 p.m.	CLUB MEETING	CARRONADE ROOM 2

AGENDA – SUNDAY, FEBRUARY 23, 2020

7:00 a.m.	SAFETY CELL (until 9:00 a.m. or later if needed.)	POOLSIDE
8:00 a.m.	OPC MEETING	BOARD ROOM
8:30 a.m.	INBOARD MEETING	CARRONADE ROOM 1
8:30 a.m.	INBOARD RACE CHAIRPERSON MEETING	CARRONADE ROOM 1
8:30 a.m.	VINTAGE	CARRONADE ROOM 1
9:00 a.m.	MODIFIED OUTBOARD MEETING	CARRONADE ROOM 2
10:00 a.m.	STOCK OUTBOARD MEETING	CARRONADE ROOM 2
11:20 a.m.	ALL MEETINGS ADJOURN	
11:30 a.m.	GENERAL MEETING	CARRONADE ROOM 1

Tom Johnston Region 6 Chairman (517) 451-9120 h (734) 344-9950 c

**Plan for Cell Testing (Dunk) 20/21 Feb 2021
Region 6 (APBA) Spring Meeting
Holiday Inn/Suites Toledo Southwest-Perrysburg, Ohio**

Dunk Test Procedures to meet Ohio Covid-19 Protocol

Masks to be worn at all times when not in pool. Pool personnel wearing scuba gear are exempt. Pool area off limits to spectators, only Region 6 Safety Team (Box 21) members (6) and max of (2) students will be permitted.

- Entry- temp will be taken and entered on waiver sign in sheet
- Air regulator will be given to student from sanitization station, for use throughout training session
- Mask to be placed on table by pool entry, prior to entering pool for prelim session
- Student will plug regulator into air system once in pool
- Instructor and student will maintain separation as much as is possible
- Student will disconnect from air system, exiting pool and placing mask on while proceeding to capsule area, student to keep mask on while receiving briefing from instructor (who will be wearing M95 mask).
- Following briefing, students will place their mask on table and enter capsule. Once seated they will plug their regulator into capsule air system.
- Following capsule egress, student will exit pool, retrieve their mask and proceed to check in area to return regulator for sanitization.
- Student is given capsule training certificate and exits pool area.

Regulators will be submersed in approved disinfectant solution after each use. Tables will be sanitized/wiped after each use.

Pool is chlorinated water, that will help in the sanitizing process.

Air tanks are 3000k psi with step down to 150psi on demand scuba regulators.